

81020303

Bipolar Forceps with Irrigation, bayonet



Technical specifications

Dimensions and Weights

Total length	200 mm
Total length in inches	7,87 "
Tip width	1,5 mm

Basic Specifications

Forceps shape	bayonet
Tip shape	straight
Connector	Flat Plug
with rinse	Yes

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418707350