


85717151

BiTech Bipolar Scissor



Technical specifications

Dimensions and Weights

Total length	175 mm
Total length in cm	17,5 cm

Basic Specifications

Scissor blade	Micro
Scissor handle	curved

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 300 Vp
Packaging unit	1
GTIN 14	04250418733342