

85620420

MICRODOME Bipolar Coagulation Electrode



Technical specifications

Dimensions and Weights

Total length	130 mm
Total length in cm	13 cm
Diameter	2 mm

Basic Specifications

Electrode shape	Dome
Ausführung	angled 20°

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 250 Vp
Packaging unit	1
GTIN 14	04250418736596