


**82460204**

**POWERGRIP 3.0 micro Bipolar Coagulation Forceps, Shaft**



## Technical specifications

### Dimensions and Weights

Total length	200 mm
Diameter	3 mm

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

### Additional Product Information

Packaging unit	1
GTIN 14	04250418742504